

H I S T O C O R E

# Configuration Guide

Original Validated HistoCore Microtome  
Accessories for Research and Industry  
Sectioning Applications




For Research Use Only. Not for use in diagnostic procedures.

**Leica**  
BIO SYSTEMS


## YOU WANT TO SECTION NON-HUMAN BIOLOGICAL SPECIMENS EMBEDDED IN PARAFFIN, UNDECALCIFIED BONE, AUTOMOTIVE OR PLASTIC PARTS?

Our reliable and precise microtomes, in combination with the original validated accessories, are designed to provide high quality sections and enhanced efficiency, while maintaining a safe and healthy workplace.

Choose the rotary microtome that lets you get the best possible section from every block and broaden your research with a wide selection of blade and specimen holders, enabling you to discover new breakthroughs in research for biomedical to industrial applications.




Welding Seam



Distribution of Active Particles in Absorber Material



Clear Coat Defect



Potential Rupture Point



TECHNICAL DATA	HistoCore BIOCUT R	HistoCore MULTICUT R	HistoCore AUTOCUT R	HistoCore NANOCUT R
Type of microtome	Mechanical	Semi Motorized	Fully Motorized	Fully Motorized
<b>GENERAL</b>				
Nominal supply voltages:	N/A	100/120/230/240 V AC	100/120/230/240 V AC	100/120/230/240 V AC
Nominal frequency:	N/A	50/60 Hz	50/60 Hz	50/60 Hz
<b>DIMENSIONS AND WEIGHT</b>				
Width (including handwheel and coarse feed wheel) Depth (including the section waste tray) Height (without top tray) W x D x H:	477 mm x 620 mm x 295 mm	477 mm x 620 mm x 295 mm	477 mm x 620 mm x 295 mm	415 mm x 620 mm x 295 mm
Weight (without accessories):	Approx. 31 kg	Approx. 31 kg	Approx. 40 kg	Approx. 40 kg
<b>MICROTOME</b>				
Section thickness setting range:	1 - 60 µm	0.5 - 100 µm	0.5 - 100 µm	0.25 - 50 µm
Trimming section thickness setting range:	10 µm, 30 µm	1 - 600 µm	1 - 600 µm	1 - 300 µm
Specimen feed:	Approx. 24 mm ±2 mm	Approx. 24 mm ±1 mm	Approx. 24 mm ±1 mm	Approx. 24 mm ±1 mm
Vertical stroke:	70 mm ±1 mm	70 mm ±1 mm	70 mm ±1 mm	70 mm ±1 mm
Maximum specimen size (H x W x D):	Large standard clamp: 55 x 50 x 30 mm Super Cassette clamp: 68 x 48 x 15 mm	Large standard clamp: 55 x 50 x 30 mm Super Cassette clamp: 68 x 48 x 15 mm	Large standard clamp: 55 x 50 x 30 mm Super Cassette clamp: 68 x 48 x 15 mm	Large Standard Clamp: 55 x 50 x 30 mm Super Cassette clamp: 68 x 48 x 15 mm
Unique force balance system	Yes	Yes	Yes	Yes
Specimen retraction:	Approx. 40 µm; can be turned off	5 - 100 µm in 5 µm increments; can be turned off	5 - 100 µm in 5 µm increments; can be turned off	5 - 50 µm (in 5 µm increments); can be turned off
<b>COARSE FEED AND MOTORIZED SECTIONING SPEEDS</b>				
Slow forward and backward speed Fast forward speed Fast backward speed (fast homing)	N/A	300 µm/s 800 µm/s 1800 µm/s	300 µm/s 800 µm/s 1800 µm/s	150 µm/s 400 µm/s 900 µm/s
Sectioning speed:	N/A (manual)	N/A (manual)	0 - 420 mm/s ±10 %	0 - 195 mm/s ±10%
Personalized coarse feed wheel	User selectable	User selectable	User selectable	N/A
Specimen orientation with zero position horizontal / vertical rotation:	± 8° / ± 8°	± 8° / ± 8°	± 8° / ± 8°	± 8° / ± 8°
Waste tray	Standard	Standard	Standard	Standard

For Research Use Only. Not for use in diagnostic procedures.

# CONFIGURE YOUR OWN RESEARCH MICROTOME


## SELECT YOUR MICROTOME

BIOCUT R	MULTICUT R	AUTOCUT R	NANOCUT R
Mechanical, manual 14 0521 58201	Semi-motorized, manual 14 0522 58221	Motorized 14 0523 58241	Motorized 14 0524 58261

## ORIENTING OR NON-ORIENTING FIXTURE FOR SPECIMEN CLAMPS

			
Rigid 14 0502 38160	Fine Directional 14 0502 37717	Directional 14 0502 38949	Segment Arc (direct assembly) 14 0502 40314


## SPECIMEN HOLDERS AND CLAMPS




- BIOCUT R only
- MULTICUT R only
- AUTOCUT R only
- NANOCUT R only or recommended



## KNIFE AND BLADE HOLDERS



## OTHER ACCESSORIES



## HOW TO CONFIGURE YOUR MICROTOME:


- First select the microtome type you are looking for: You have the choice between 2 manual and 2 fully automated microtomes.
- Decide whether you would like to use a specimen orientation (fine directional, directional or via segment arc) or select a fixed rigid object head for extra stability.
- Are you sectioning cassettes (standard or super size), square or round specimen blocks, or EM (Semi-thin for Electron Microscope) specimens? Select the specimen holder which is ideal for your specimen.
- Depending on the type of blade or knife you will use, select the blade or knife holder.
- The knife holder base and Knife holder GD for glass and diamond knives can additionally be equipped with a backlighting system, which helps you to better align the specimen to the knife.
- The color codes (■ ■ ■ ■) indicate with which instrument the accessory can be used. Accessories without color codes can be used with all 4 microtomes.
- A large variety of disposable blades or reusable knives can be found in the consumable section of the [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com) Website.

# NON-HUMAN PARAFFIN SECTIONING APPLICATIONS


## SELECT YOUR MICROTOME




## ORIENTING OR NON-ORIENTING FIXTURE FOR SPECIMEN CLAMPS



## CASSETTES OR BLOCKS



## DISPOSABLE BLADES OR REUSABLE KNIVES



## OTHER ACCESSORIES



- BIOCUT R only
- MULTICUT R only
- AUTOCUT R only
- NANOCUT R only

- Do you want to work manually or with an automated microtome? Select Your Microtome.
- To have both hands free for section collection, the AUTOCUT R can be equipped with a footswitch to easily start and stop the motorized sectioning.
- Is a specimen orientation required (e.g. for re-cuts) and how precise does it need to be? You have the choice of either the fine directional fixture for specimen clamps with zero indicators and click stops every 2° or the directional fixture for specimen clamps, with 8° XY orientation. Both orientation systems can be equipped with the quick clamping system, allowing rapid specimen clamp exchange or clamp removal for cleaning.
- If you prefer a fixed and stable object head, the rigid fixture (with the integrated quick clamping system) is the right choice for you.
- Are you going to section specimens embedded in standard or super size cassettes or paraffin blocks? Do you want to cool your standard cassettes

while sectioning to keep the block colder longer, especially when sectioning serial sections? Find the appropriate clamp for your specimen in the cassettes or blocks section.

- Depending on whether you prefer disposable blades or reusable knives, select the blade or knife holder. The 2-in-1 Blade Holder (to be used with the blade holder base) can be used with either high or low profile blades.
- Knife holder E with a trough is used with the knife holder base and low profile blades. The trough can be filled with water and the sections will stretch on the water surface and can be collected with an object slide.
- When using reusable c or d profile knives, Knife Holder N is the best choice.
- The Knife Holder Base (for Knife Holder N and E with a trough) can be equipped with a backlighting system, which helps you to better align the specimen to the knife.
- A large variety of disposable blades or reusable knives can be found in the consumable section of the [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com) Website.



Knife Holder E with Trough

#	Recommendation	Order No.
	BIOCUT R Basic Instrument	14 0521 58201
OR	MULTICUT R Basic Instrument	14 0522 58221
OR	AUTOCUT R Basic Instrument	14 0523 58241
1	Fine-Directional Fixture for Specimen Clamps	14 0502 37717
2	Quick Clamping System	14 0502 37718
3	Leica RM CoolClamp	14 0502 46573
4	Blade Holder Base	14 0502 55546
5	2-in-1 Blade Holder	14 0502 54497
6	Low profile disposable blades e.g.	14 0358 38382
	High profile disposable blades e.g.	14 0358 38383



# CONFIGURE YOUR OWN MICROTOME FOR SEMI-THIN SECTIONING AND 3D APPLICATIONS

SELECT YOUR MICROTOME

ORIENTING OR NON-ORIENTING FIXTURE FOR SPECIMEN CLAMPS

APPLICATION

## NANOCUT R

Motorized  
14 0524 58261

Rigid  
14 0502 38160

Fine Directional  
14 0502 37717

Directional  
14 0502 38949

Segment Arc (direct assembly)  
14 0502 40314

Quick Clamping System  
14 0502 37718

Segment Arc (indirect assembly)  
14 0502 29969

3D

SEMI-THIN SECTIONING

Round Specimen Holder with 3 inserts: 6, 15 and 25 mm ø  
14 0502 38002

EM Specimen Holder Fixture  
14 0502 29968

EM Flat Specimen Holder  
14 0355 10405

EM Universal Specimen Holder  
14 0356 10868

Knife Holder Base  
14 0502 37962

Backlighting  
14 0502 38719

Special Wrench  
14 0356 10869

Knife Holder E-TC  
14 0502 37997

Knife Holder GD  
14 0502 39052

Disposable Tungsten Carbide Blades TC-65  
14 0216 26379

Diamond Knives  
Triangular Glass Knives

OTHER ACCESSORIES


Foot Switch  
14 0502 38257

Microscope Carrier with LED Hi Power Spots and Microscope  
14 0502 38463

Combined Carrier with Magnifier and Illumination  
14 0502 29971 120V/60Hz  
14 0502 29972 230V/50Hz

LED 1000 Hi-Power Spots, 2 Arms  
14 6000 04826

LED 1000 Control Unit  
14 6000 04825



Raw data from a E15,5 stage mouse embryo



### 3D APPLICATION – FOLLOW THE MERLOT DECISION TREE

- The NANOCUT R is the instrument of choice for 3D reconstruction with its precise stop in the object head's optimal position for photo-capture of the cut sample surface when running 3-dimensional reconstruction programs.
- Choose whether you would like to use a specimen orientation (fine directional or directional) or prefer a fixed rigid object head for extra stability.
- The quick clamping system will allow you to quickly remove the specimen clamp.
- As most specimens are embedded using molds with round specimen clamp adapters, the use of a round specimen clamp is recommended.
- For 3-dimensional reconstruction programs, high quality surfaces of the cut block surface are required. This can be obtained by using disposable tungsten carbide TC-65 blades, which are clamped in the E-TC blade holder.

### SEMI-THIN SECTIONING – FOLLOW THE GREEN DECISION TREE


- The NANOCUT R is the instrument of choice for semi-thin sectioning with a minimum section thickness of 250 nm. The use of a footswitch is recommended allowing you to have both hands free for the section preparation.
- The specimen size for semi-thin sectioning is limited. Specimens are embedded in resin and will be clamped in either an EM flat or universal specimen holder. The specimens are clamped in the holders using the special wrench.
- Often a segment arc is used to clamp the EM holders – you have the choice to either select the segment arc, which is directly clamped onto the cylinder of the microtome, which gives you the most stability, or the segment arc with indirect assembly, which is clamped on either the fine directional or directional fixture for specimen clamps.
- Semi-thin sections can be produced using either glass or diamond knives, which are clamped in the Knife holder GD. Backlighting is recommended to help better align the specimen to the knife.
- LED illumination and the microscope will help you to prepare and collect the section ribbons.

### RECOMMENDED SOLUTION FOR 3D RECONSTRUCTION

#	Recommendation	Order No.
	NANOCUT R Basic Instrument	14 0521 58261
1	Fine-Directional Fixture for Specimen Clamps	14 0502 37717
2	Quick Clamping System	14 0502 37718
3	Round specimen holder	14 0502 38002
4	Knife holder base	14 0502 37962
5	Backlighting	14 0502 38719
6	Knife holder E-TC for TC-65 disposable blades	14 0502 37997
7	TC-65 disposable blades	14 0216 26379
8	Foot switch	14 0502 38257

### RECOMMENDED SOLUTION FOR SEMI-THIN SECTIONING

#	Recommendation	Order No.
	NANOCUT R Basic Instrument	14 0521 58261
1	Fine-Directional Fixture for Specimen Clamps	14 0502 37717
2	Quick Clamping System	14 0502 37718
3	EM specimen holder fixture	14 0502 29968
4	EM holder flat or	14 0355 10405
5	EM holder universal	14 0356 10868
6	Special wrench	14 0356 10869
7	Knife holder GD	14 0502 39052
8	Diamond knife	N/A
9	Microscope carrier w/LED Hi-Power spots	14 0502 38463
10	LED 1000 control unit	14 6000 04825
11	Foot switch	14 0502 38257



# CONFIGURE YOUR OWN MICROTOME FOR SPECIMENS OF DIFFERENT HARDNESS AND DESIRED SECTION THICKNESS OR SURFACE QUALITY

## SELECT YOUR MICROTOME


## ORIENTING OR NON-ORIENTING FIXTURE FOR SPECIMEN CLAMPS

## SPECIMEN SIZE AND FORMAT

## SPECIMEN HARDNESS SECTION OR SURFACE

## DESIRED MINIMAL SECTION THICKNESS

## OTHER ACCESSORIES



To section resin embedded biological or plastic specimens of different hardness, an automated microtome is recommended.

- Select either the AUTOCUT R for larger blocks and section thicknesses down to 0.5 µm or the NANOCUT R for specimens requiring extreme slow cutting speeds, thinner sections down to 250 nm or high-quality surfaces of specimens.
- To have both hands free for section collection, select a footswitch to easily start and stop motorized sectioning.
- Is a specimen orientation required (e.g. for target preparation) and how precise does it need to be?
- The fine directional fixture for specimen clamps with zero indicators and click stops every 2°.
- The directional fixture for specimen clamps, with 8° XY orientation, with the opportunity to turn the clamp by 360° to reduce cutting force.
- Both orientation systems are equipped with the quick clamping system, for rapid clamp exchange.
- If you prefer a fixed and stable object head, the rigid fixture (with the integrated quick clamping system) is the right choice for you.
- In order to prevent chemical or thermal effects on the sample to be cut, many users do not embed the plastic sample in synthetic resin or paraffin and prefer direct clamping.
- For square or angular specimens one of the standard specimen clamps is recommended.

## RECOMMENDED SOLUTION FOR RESIN SECTIONING E.G. BONE OR HARDER INDUSTRIAL SAMPLES

#	Recommendation	Order No.
	AUTOCUT R Basic Instrument	14 0523 58241
1	Fine-Directional Fixture for Specimen Clamps	14 0502 37717
2	Quick Clamping System	14 0502 37718
3	Round specimen holder or	14 0502 38002
4	Standard clamp 40 x 40 mm	14 0502 37998
5	Knife holder base	14 0502 37962
6	Backlighting	14 0502 38719
7	Knife holder NZ for reusable tungsten carbide (TC) knives	14 0502 39052
8	16 cm TC knife d profile	14 0216 04813
or 9	Knife holder E-TC for disposable TC blades	14 0502 37997
10	TC-65 disposable blades	14 0216 26379
11	Foot switch	14 0502 38257




- To stabilize foils while sectioning, it is mounted like a sandwich in between two plastic parts. The sandwich is clamped in the foil clamp, which is mounted into the standard clamp.
- For round specimens of 5, 18 or 25 mm diameter, the round specimen clamp can be used. For other diameters use one of the standard clamps with the V-insert.
- For smaller flat or capsule size specimens, choose the EM flat or universal specimen holder. They are clamped into the EM holder fixture.
- What is the hardness of your specimen? Do you want to produce sections or are you going to investigate the cut surface of the remaining block? Which minimum section thickness are you looking for?
- For soft plastics, such as polyethylene and a desired section thickness of approx. 5 µm, use disposable blades and the 2-in-1 blade holder.
- For plastics, such as Glycol methacrylate (GMA e.g. HistoResin) or hard plastics, such as Methyl methacrylate (MMA) or Epoxy resin, the use of a tungsten carbide blade/knife is suitable.
- Down to 2 µm: Disposable TC-65 blade mounted in the E-TC knife holder.
- For thinner sections, a reusable tungsten carbide knife (with d profile) clamped in the NZ knife holder should be used.
- For high quality surfaces, e.g for defect analysis, either the tungsten carbide blade/knife can be used, however, for highest possible surface quality, a Diamond knife with GD knife holder, is recommended.
- Magnifier or microscope with LED illumination will help you to prepare and collect the sections.

## RECOMMENDED SOLUTION FOR SOFT PLASTIC PART SECTIONING (INDUSTRY)


#	Recommendation	Order No.
	AUTOCUT R Basic Instrument	14 0523 58241
1	Fine-Directional Fixture for Specimen Clamps or	14 0502 37717
2	Directional fixture for Specimen Clamps	14 0502 38949
3	Quick Clamping System	14 0502 37718
4	Standard clamp 50 x 55 mm	14 0502 38005
5	Foil Clamp	14 0402 09307
6	V-Insert	14 0502 38000
7	Blade Holder Base	14 0502 55546
8	2-in-1 Blade Holder	14 0502 54497
9	Low profile disposable blades e.g. 819	14 0358 38382
	High profile disposable blades e.g. 819	14 0358 38383
10	Foot switch	14 0502 38257




# RESEARCH MICROTOMES




Impurities in a foil 50 µm




Clearcoat inhomogeneities 20 µm




Sandwich foil 5 µm



Proximal tibia of a rat 5 µm



Huge glass fiber perforates the topcoat 20 µm



Plastic form 20 µm

## CONFIGURED HISTOCORE ROTARY MICROTOMES FOR YOUR NON-HUMAN PARAFFIN SECTIONING RESEARCH APPLICATION

	ORDER NUMBER	ORDER NUMBER	ORDER NUMBER
Standard delivery includes	149BIOR00C1	149MULTIRC1	149AUTOR0C1
Precision orientation: Fine directional fixture for specimen clamps	✓	✓	✓
Quick clamping system	✓	✓	✓
Universal cassette clamp	✓	✓	✓
Blade holder base	✓	✓	✓
2-in-1 blade holder	✓	✓	✓
Standard waste tray	✓	✓	✓
Top tray	✓	✓	✓

LeicaBiosystems.com

Copyright © 2019 Leica Biosystems Imaging, Inc. All Rights Reserved. LEICA and the Leica logo are registered trademarks of Leica Microsystems IR GmbH. Aperio is a trademark of the Leica Biosystems group of companies in the USA and optionally in other countries. Other logos, product and/or company names might be trademarks of their respective owners.

190671 Rev A 08/2019